

# Roc


Roc uses lateral movement to control a the vehicle to avoid oncoming traffic. It uses 2 markers to drive the car left or right to see how long the player can make it.



Step 1:  
Download and install app.  
Follow onscreen instructions and request license from [license@idoneusdigital.com](mailto:license@idoneusdigital.com).



Step 2:  
Download and install IdoneusConnect.  
Follow onscreen instructions and request license from [license@idoneusdigital.com](mailto:license@idoneusdigital.com)  
Place all licenses in C:\user\xxxx\IdoneusDigital\



Step 3:  
Attached Required markers (Left Hip & Right Hip) and begin streaming live marker labels.  
(Not required for use with MS Kinect)



Step 4:  
Start IdoneusConnect and confirm system is streaming frames. (see image)

# Roc



## Step 5:

Start App and confirm it is connected by verifying the small “Plug in” icon in the bottom left corner is in the “Plugged In” state (can take a few seconds). Then hit the gear “Settings” icon in the bottom right corner to bring up the core settings screen.



## Step 6:

Retarget Motion Capture markers to App markers. Roc requires 2 labeled markers to operate correctly however you can use as many as needed in order to get marker naming to stream properly. (Clusters)

Once connected set the correct orientation of the game so it projects correctly onto wall with the projector or TV



## Step 7:

Purchase/choose vehicle and select environment to race.



## Step 8:

The longer the player goes the more money they collect to unlock better vehicles.