PRINCIPLES OF INSTRUMENTAL ANALYSIS

FIFTH EDITION

SKOOG · HOLLER · NIEMAN

PRINCIPLES OF Instrumental Analysis

Fifth Edition

Douglas A. Skoog

Stanford University

F. James Holler

University of Kentucky

Timothy A. Nieman

University of Illinois at Urbana-Champaign

#
THOMSON LEARNING

Australia • Canada • Mexico • Singapore • Spain United Kingdom • United States

BROOKS/COLE THOMSON LEARNING

Chemistry Editor: John Vondeling

Developmental Editors: Marc Sherman, Jennifer Bortel

Project Editor: Frank Messina

Marketing Strategist: Angus McDonald

Production Manager: Charlene Catlett Squibb

Text Designer: Ruth Hoover Art Director: Lisa Adamitis Copy Editor: Vince Pecuilis Cover Designer: Chazz Bjanes

Cover Image: John T. Foley, Department of Physics and Astronomy, Mississippi State University

Cover Printer: Lehigh Press, Inc. Compositor: The Clarinda Company

Printer: R.R. Donnelley—Crawfordsville, MD

Copyright © 1998, 1992, 1985, 1980, 1971 Thomson Learning, Inc. Thomson Learning™ is a trademark used herein under license.

ALL RIGHTS RESERVED. No part of this work covered by the copyright hereon may be reproduced or used in any form or by any means—graphic, electronic, or mechanical, including but not limited to photocopying, recording, taping, Web distribution, information networks, or information storage and retrieval systems—without the written permission of the publisher.

Printed in the United States of America 10 11 12 05 04 03 02

For more information about our products, contact us at:

Thomson Learning Academic Resource Center 1-800-423-0563

For permission to use material from this text, contact us by:

Phone: 1-800-730-2214 **Fax:** 1-800-730-2215 **Web:** http://www.thomsonrights.com

Library of Congress Catalog Card Number: 97-67376

PRINCIPLES OF INSTRUMENTAL ANALYSIS, Fifth Edition

ISBN: 0-03-002078-6

Asia

Thomson Learning 60 Albert Street, #15-01 Albert Complex Singapore 189969

Australia

Nelson Thomson Learning 102 Dodds Street South Melbourne, Victoria 3205 Australia

Canada

Nelson Thomson Learning 1120 Birchmount Road Toronto, Ontario M1K 5G4 Canada

Europe/Middle East/Africa

Thomson Learning Berkshire House 168-173 High Holborn London WC1 V7AA United Kingdom

Latin America

Thomson Learning Seneca, 53 Colonia Polanco 11560 Mexico D.F. Mexico

Spain

Paraninfo Thomson Learning Calle/Magallanes, 25 28015 Madrid, Spain

Standards

he two ws the lecular acking a chro-S G n on a nber of ize-exof glujuices. 00, was ophilic ng con-G n chroolecular er polyimes of or a see same

Cranberry

Orange

Figure 28-29 Determination of glucose (G), fructose (F), and sucrose (S) in canned juices. (Courtesy of Perkin-Elmer Corporation, Norwalk, CT.)

Time, min

Apple

Pineapple

The disadvantages are: (1) only a limited number of bands can be accommodated because the time scale of the chromatogram is short; and (2) inapplicability to samples of similar size, such as isomers. Generally, a 10% difference in molecular weight is required for reasonable resolution.

clusion separa-

 V_o and

narrow

m from the stativation cking.

28H THIN-LAYER CHROMATOGRAPHY

Planar chromatographic methods include thin-layer chromatography (TLC) and paper chromatography (PC), and electrochromatography. Each makes use of a flat, relatively thin layer of material that is either selfsupporting or is coated on a glass, plastic, or metal surface. The mobile phase moves through the stationary phase by capillary action, sometimes assisted by gravity or an electrical potential. Planar chromatography is sometimes called two-dimensional chromatography, although this description is not strictly correct inasmuch as the stationary phase does have a finite thickness.

Currently, most planar chromatography is based upon the thin-layer technique, which is faster, has bet-

ter resolution, and is more sensitive than its paper counterpart. This section is devoted to thin-layer methods only.

28H-1 The Scope of Thin-Layer Chromatography

In terms of theory, the types of stationary and mobile phases, and applications, thin-layer and column liquid chromatography are remarkably similar. In fact, an important use of thin-layer chromatography, and the reason this topic is included here, is to serve as a guide to the development of optimal conditions for performing separations by column liquid chromatography. The advantages of following this procedure are the speed and low cost of the exploratory thin-layer experiments. Some chromatographers, in fact, take the position that thin-layer experiments should always precede column experiments.

In addition to its use in developing column chromatographic methods, thin-layer chromatography has become the workhorse of the drug industry for the all-important determination of product purity. It has also found widespread use in clinical laboratories and is the backbone of many biochemical and biological studies. Finally, it finds widespread use in the industrial laboratories. Al As a consequence of these many areas of application, it has been estimated that at least as many analyses are performed by thin-layer chromatography as by high-performance column liquid chromatography. Al

28H-2 How Thin-Layer Separations Are Performed

Typical thin-layer separations are performed on flat glass or plastic plates that are coated with a thin and adherent layer of finely divided particles; this layer constitutes the stationary phase. The particles are similar to those described in the discussion of adsorption, normal-and reversed-phase partition, ion-exchange, and size-exclusion column chromatography. Mobile phases are also similar to those employed in high-performance column liquid chromatography.

Thin-Layer Plates

Thin-layer plates are available from several commercial sources at costs that range from \$1 to \$10 per plate. The common plate sizes in centimeters are $5\times 20,\,10\times 20,\,$ and $20\times 20.\,$ Commercial plates come in two categories, conventional and high-performance. The former have thicker layers (200 to 250 μm) of particles having nominal particle sizes of 20 μm or greater. High-performance plates usually have film thicknesses of 100 μm and particle diameters of 5 μm or less. High-performance plates, as their name implies, provide sharper separations in shorter times. Thus, a conventional plate typically will exhibit 2000 theoretical plates in 12 cm with a development time of 25 min. The corresponding figures for a high-performance plate are 4000 theoreti-

cal plates in 3 cm requiring 10 min for development. High-performance plates suffer from the disadvantage of having a significantly smaller sample capacity.


Sample Application

Sample application is perhaps the most critical aspect of thin-layer chromatography, particularly for quantitative measurements. Usually the sample, as a 0.01% to 0.1% solution, is applied as a spot 1 to 2 cm from the edge of the plate. For best separation efficiency, the spot should have a minimal diameter—about 5 mm for qualitative work and smaller for quantitative analysis. For dilute solutions, three or four repetitive applications are used, with drying in between.

Manual application of samples is performed by touching a capillary tube containing the sample to the plate or by use of a hypodermic syringe. A number of mechanical dispensers, which increase the precision and accuracy of sample application, are now offered commercially.

Plate Development

Plate development is the process in which a sample is carried through the stationary phase by a mobile phase; it is analogous to elution in liquid chromatography. The most common way of developing a plate is to place a drop of the sample near one edge of the plate and mark its position with a pencil. After the sample solvent has evaporated, the plate is placed in a closed container saturated with vapors of the developing solvent. One end of the plate is immersed in the developing solvent, with care being taken to avoid direct contact between the sample and the developer (Figure 28-30). The developing solvent travels up the plate, being drawn by capillary action between the fine particles. As the developing solvent travels past the point of sample application, it dissolves



⁴² T. H. Mauch II, Science, 1982, 216, 161; J. C. Touchstone, LC-GC, 1993, 11(6), 404.




Figure 28-30 A typical development chamber.