

COMPLEX-AID is a source of nitrogen for use in agronomic crops.

COMPLEX AID™ 2-0-0

GUARANTEED ANALYSIS

Total Nitrogen (N) 2.0%

Derived from Hydrolyzed Brewer's Yeast

STORAGE & PRECAUTIONS: Keep container tightly sealed when not in use. Store in a cool, dry place between 34° and 100° F. Do not store in direct sunlight. Do not store diluted product. Liquid should not be consumed. When used in conjunction with fungicides and insecticides, the following precautions apply: Do not inhale mist. Wash with soap and water after use or skin contact. Prevent contamination of food, utensils and drinking water. Dispose of contents/container to hazardous or special waste collection point in accordance with local, regional, national and/or international regulation. **KEEP OUT OF REACH OF CHILDREN. SHAKE WELL BEFORE USE.**

Specimen Label

Alltech®
CROP SCIENCE

Improcrop U.S.A. Inc. d/b/a
Alltech Crop Science
3031 Catnip Hill Road,
Nicholasville, KY 40356
Tel: (859) 885-9613
REV 04.16 SMB

2.5 Gal.
Net weight 26.9 lbs (12.2 kg)

Manufacture date: XXXXXXXX
Batch number: XXXXXXXX