						💄 Sign i	n Create an account
Europe PMC	About	Tools	Developers	Help	Europe PMC plus		
Search worldwide, life-sciences literature							
CITES:17960406_MED						Q Search	Advanced Search
Coronavirus articles and preprints Searc	h examples:	"breast can	cer" <u>Smith J</u>				
Recent history Saved searches							

Fluorescence characteristics and pharmacokinetic properties of a novel self-adhesive 5-ALA patch for photodynamic therapy of actinic keratoses.

Fauteck JD¹, Ackermann G, Birkel M, Breuer M, Moor AC, Ebeling A, Ortland C

Author information >

Archives of Dermatological Research, 25 Oct 2007, 300(2) 53-60 DOI: 10.1007/s00403-007-0787-0 PMID: 17960406

Share this article 🛛 😏 🛅 🧍

Abstract

Actinic keratosis (AK) can be treated by photodynamic therapy (PDT), which is becoming a well-established tool in dermatology. Normally a precursor of the photosensitiser is applied topically and converted into protoporphyrin IX (PPIX) in the cells. By activating PPIX with light, the dysplastic cells will be destroyed. We report the results of two clinical studies investigating the properties of a novel self-adhesive 5-ALA-patch (PD P 506 A) intended for PDT of m to moderate AK on the face and head. The studies investigated the influence of patch application duration on PPIXspecific fluorescence and the pharmacokinetic properties of the 5-ALA patch. The PPIX fluorescence in AK lesions a normal skin after patch application (intraindividual comparison; application for 2, 3, 4, 5 h) was investigated in 13 patients using DYADERM Professional (Biocam). In the subsequent pharmacokinetic study 12 patients were treated with 8 patches each (4 h application). 5-ALA and PPIX were analysed in plasma (over 24 h) and urine (over 12 h). PPI specific fluorescence measured immediately after patch removal increased with increasing application duration to maximum at 4-h application. The fluorescence in AK lesions was more intense than in normal skin. A small increase of 5-ALA plasma concentrations was observed in 10 of 12 patients after applying 8 patches for 4 h, which rapidly declined to normal values after patch removal. The maximum increase was 3.7-fold of the pre-dose 5-ALA plasma concentration. No PPIX-concentrations above the lower limit of quantification were observed. PPIX-specific fluorescence in AK lesions can be steered by application duration of this novel 5-ALA patch. Application is safe and well tolerable. The observed small rise in 5-ALA plasma concentrations is regarded clinically irrelevant. Clinical efficacy of the patch in PDT will be investigated in further clinical trials.

Full text links

References

Citations & impact

Impact metrics

This website requires cookies, and the limited processing of your person. By using the site you are agreeing to this as outlined in our <u>privacy notice</u> Biofrontera Exhibit 1034 Biofrontera Inc. et al. v. DUSA Pharmaceuticals, Inc. IPR2022-00056

						💄 Sign i	n Create an account
Europe PMC	About	Tools	Developers	Help	Europe PMC plus		
Search worldwide, life-sciences literature							
CITES:17960406_MED						Q Search	Advanced Search
Coronavirus articles and preprints Searc	h examples:	"breast can	cer" <u>Smith J</u>				
Recent history Saved searches							

Fluorescence characteristics and pharmacokinetic properties of a novel self-adhesive 5-ALA patch for photodynamic therapy of actinic keratoses.

Fauteck JD¹, Ackermann G, Birkel M, Breuer M, Moor AC, Ebeling A, Ortland C

Author information >

Archives of Dermatological Research, 25 Oct 2007, 300(2) 53-60 DOI: 10.1007/s00403-007-0787-0 PMID: 17960406

Share this article 🛛 🔰 🛅 🛉

Abstract

Actinic keratosis (AK) can be treated by photodynamic therapy (PDT), which is becoming a well-established tool in dermatology. Normally a precursor of the photosensitiser is applied topically and converted into protoporphyrin IX (PPIX) in the cells. By activating PPIX with light, the dysplastic cells will be destroyed. We report the results of two clinical studies investigating the properties of a novel self-adhesive 5-ALA-patch (PD P 506 A) intended for PDT of m to moderate AK on the face and head. The studies investigated the influence of patch application duration on PPIXspecific fluorescence and the pharmacokinetic properties of the 5-ALA patch. The PPIX fluorescence in AK lesions a normal skin after patch application (intraindividual comparison; application for 2, 3, 4, 5 h) was investigated in 13 patients using DYADERM Professional (Biocam). In the subsequent pharmacokinetic study 12 patients were treated with 8 patches each (4 h application). 5-ALA and PPIX were analysed in plasma (over 24 h) and urine (over 12 h). PPI specific fluorescence measured immediately after patch removal increased with increasing application duration to maximum at 4-h application. The fluorescence in AK lesions was more intense than in normal skin. A small increase of 5-ALA plasma concentrations was observed in 10 of 12 patients after applying 8 patches for 4 h, which rapidly declined to normal values after patch removal. The maximum increase was 3.7-fold of the pre-dose 5-ALA plasma concentration. No PPIX-concentrations above the lower limit of quantification were observed. PPIX-specific fluorescence in AK lesions can be steered by application duration of this novel 5-ALA patch. Application is safe and well tolerable. The observed small rise in 5-ALA plasma concentrations is regarded clinically irrelevant. Clinical efficacy of the patch in PDT will be investigated in further clinical trials.

Full text links

References

Citations & impact

Impact metrics

This website requires cookies, and the limited processing of your personal data in order to function. By using the site you are agreeing to this as outlined in our <u>privacy notice</u> and <u>cookie policy</u>.



Mentioning Ø 27 Contrasting ⑦ 0

Article citations

Smartphone fluorescence imager for quantitative dosimetry of protoporphyrin-IX-based photodynamic therapy in skin. Ruiz AJ, LaRochelle EPM, Gunn JR, Hull SM, Hasan T, Chapman MS, Pogue BW J Biomed Opt, 25(6):1-13, 01 Dec 2019 Cited by: 1 article | PMID: 31820594 | PMCID: PMC6901011 Free to read & use

Photodynamic therapy simplified: nonprepared, moderate-grade actinic keratosis lesions respond equally well to 5-aminolaevulinic acid patch photodynamic therapy as do mild lesions. Szeimies RM, Hauschild A, Ortland C, Moor AC, Stocker M, Surber C Br J Dermatol, 173(5):1277-1279, 26 Aug 2015 Cited by: 3 articles | PMID: 25940103

PpIX fluorescence combined with auto-fluorescence is more accurate than PpIX fluorescence alone in fluorescence detection of non melanoma skin cancer an intra patient direct comparison study van der Beek N de Leeuw J Demmendal C Bjerring P Neumann HA Lasers Surg Med, 44(4):271-276, 13 Dec 2011 Cited by: 6 articles | PMID: 22170313

New aspects in photodynamic therapy of actinic keratoses. Smits T, Moor AC J Photochem Photobiol B 96(3) 159 169 13 Jun 2009 Cited by 5 articles | PMID 19592269 Review

An overview of topical photodynamic therapy in dermatology. Ibbotson SH Photodiagnosis Photodyn Ther, 7(1):16-23, 29 Dec 2009 Cited by 17 articles LPMID 20230989

This website requires cookies, and the limited processing of your personal data in order to function. By using the site you are agreeing to this as outlined in our <u>privacy notice</u> and <u>cookie policy</u>.

I agree, dismiss this banner

Abstract Full text 🗹 References

Citations & impact Similar Articles

002

Similar Articles

About

About Europe PMC Funders Joining Europe PMC Governance Roadmap Outreach

Tools

Tools overview ORCID art cle cla m ng Journal list Grant finder External links service RSS feeds Annotations Annotat ons subm ss on serv ce

Developers

Developer resources Art cles RESTful API Grants RESTful API API case stud es SOAP web service Annotat ons API OAI service Bulk downloads Developers Forum

Help

Help using Europe PMC Search synta reference Contact us C He Fe Ti B Te

D

Let u know how we are doing



Europe PMC is part of the ELIXIR infrastructure Europe PMC is an ELIXIR Core Data Resource <u>Learn more ></u>

Europe PMC is a service of the Europe PMC Funders' Group, in partnership with the European Bioinformatics Institute; and in cooperation with the National Center for Biotechnology Information a includes content provided to the PMC International archive by participating publishers.



Contact us | Privacy | Terms of use | Copyright | Accessibility

This website requires cookies, and the limited processing of your personal data in order to function. By using the site you are agreeing to this as outlined in our <u>privacy notice</u> and <u>cookie policy</u>.

							👗 Sign	n Create an account
		About	Tools	Developers	Help	Europe PMC plus		
Search worldwide	e, life-sciences l	literature						
CITES:17960406_ME	Ð						Q Search	Advanced Search
Coronavirus articles an	d preprints Search	examples:	"breast can	cer". Smith J				
Recent history	Saved searches							
		1-9 of 9 re	sults					
		Sort by:	Date rece	ived	-			
		Smartpho	one fluore	scence imager f	or quanti	tative dosimetry of prot	toporphyrin-IX-based	photodynamic the
		in skin.						
Search only				EPM, Gunn JR, Hi -13, 01 Dec 2019	ull SM, Ha	san T, Chapman MS, Pog	gue BW	
Type				MID: 31820594	PMCID: PN	AC6901011		
Type ⑦		+ Add to e	export list					Free to read &
	IS (7)							
 Reviews (2) Preprints (0) 		-				ared, moderate-grade ac nic therapy as do mild le		s respond equally
					-	C, Stocker M, Surber C	.510115.	
Free full text 🔊				:1277-1279, 26 Au	ıg 2015			
Free to read (1)		Cited by: 3 + Add to e		PMID: 25940103				
Free to read & u	ISE (1)							
Date		PpIX fluo	rescence	combined with a	auto-fluor	escence is more accura	te than PpIX fluoresco	ence alone in
2019 (1)						in cancer: an intra-patie	ent direct comparisor	study.
2019 (I)				Leeuw J, Demme 4):271-276, 13 De		jerring P, Neumann HA		
2015 (1)			-	PMID: 22170313	02011			
Custom date range	e 🕨	+ Add to e	export list					
		<mark>An overvi</mark> Ibbotson		ical photodynan	nic thera	oy in dermatology.		
		-		odyn Ther, 7(1):16		ec 2009		
		Cited by: 1 + Add to e		PMID: 20230989				Rev
		- Add 10 0	export list					Key
		Long-tern data.	m follow-u	ıp of photodyna	mic thera	py with a self-adhesive	5-aminolaevulinic aci	d patch: 12 month
			RM, Stock	fleth E, Popp G, E	Borrosch F	, Brüning H, Dominicus F	R, Mensing H, Reinhold	U, Reich K, Moor
				C, Brunnert M, H		Ą		
		-		:410-414, 30 Jun 2 PMID: 19804593				
		+ Add to e		FINID, 19604595				
		New aspe	ects in pho	otodynamic the	apy of ac	tinic keratoses.		
		Smits T, N				2000		
				<mark>iol B</mark> , 96(3):159-16 PMID: 19592269	9, 13 Jun 2	2009		
		+ Add to e						Rev
This website re	quires cookies, ar	nd the limit	ed process	sing of your pers	onal data	in order to function.	Lagran dissuirs th	is houses

By using the site you are agreeing to this as outlined in our privacy notice and cookie policy.

Lasers Surg Med, 41(2):96-103, 01 Feb 2009 Cited by: 21 articles | PMID: 19226578 + Add to export list

Optimization of photodynamic therapy with a novel self-adhesive 5-aminolaevulinic acid patch: results of th randomized controlled phase III studies.

Hauschild A, Stockfleth E, Popp G, Borrosch F, Brüning H, Dominicus R, Mensing H, Reinhold U, Reich K, Moor A Stocker M, Ortland C, Brunnert M, Szeimies RM Br J Dermatol, 160(5):1066-1074, 16 Feb 2009 Cited by: 39 articles | PMID: 19222455 + Add to export list

Effective photodynamic therapy of actinic keratoses on the head and face with a novel, self-adhesive 5aminolaevulinic acid patch. Hauschild A, Popp G, Stockfleth E, Meyer KG, Imberger D, Mohr P, Itschert G, Kaufmann R, Neuber K, Frambach Gollnick H, Brunnert M, Stocker M, Ortland C, Karrer S Exp Dermatol, 18(2):116-121, 17 Jul 2008 Cited by: 21 articles | PMID: 18643849 + Add to export list

1-9 of 9 results

About

- About Europe PMC Funders Joining Europe PMC Governance Roadmap Outreach
- Tools Tools overview ORCID article claiming Journal list Grant finder External links service RSS feeds Annotations Annotations submission service
- Developers Developer resources Articles RESTful API Grants RESTful API API case studies SOAP web service Annotations API OAI service Bulk downloads Developers Forum

Help	С
Help using Europe PMC	He
Search syntax reference	Fe
Contact us	Tv
	BI
	Te
	D

Let us know how we are doing.

Europe PMC is part of the ELIXIR infrastructure Europe PMC is an ELIXIR Core Data Resource Learn more >

Europe PMC is a service of the Europe PMC Funders' Group, in partnership with the European Bioinformatics Institute; and in cooperation with the National Center for Biotechnology Information a includes content provided to the PMC International archive by participating publishers.



Contact us | Privacy | Terms of use | Copyright | Accessibility

This website requires cookies, and the limited processing of your personal data in order to function. By using the site you are agreeing to this as outlined in our privacy notice and cookie policy.